

ADHESIVES MANUFACTURERS  
WEIGHING INDIVIDUAL COMPONENTS TO DEVELOP A NEW ADHESIVE IN THE LABORATORY

PRECISION BALANCES | KERN 440

When developing a new adhesive in the laboratory, tailored to the individual needs of the customer, the individual components are precisely weighed. This is important to follow the recipe exactly and make the mix repeatable. A precision balance is ideally suited for this meticulous weighing of the components.



Model	Weighing capacity	Readability	Weighing surface (d)	Price
440-47N	2000 g	0.1 g	130×130 mm	€ 300.-
<b>Suitable Accessories:</b> <ul style="list-style-type: none"><li>• DAkkS calibration certificate 963-127 - 963-127 - €103.-</li></ul>				



<https://partner.kern-sohn.com/Segment/B040E>

\* DAkkS calibrations are also possible on site  
All prices are net, valid until 07.12.2025 plus shipping cost

Personal consultation from our product specialist:



Aleksandar Delic  
Productspecialist  
Tel : +49 7433 9933-120  
aleksandar.delic@kern-sohn.com

Application number: EN - 1036

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen, Germany  
info@kern-sohn.com  
[www.kern-sohn.com](http://www.kern-sohn.com)