

PLASTERERS AND STUCCOWORKERS
WEIGHING COLOR PIGMENTS FOR A PUTTY

PRECISION BALANCES | KERN PCB

If a tinted filler is required for the spatula technique, pigments must be stirred in. These need to be weighed carefully as the same shade needs to be achieved in a new batch. This level of accuracy can only be achieved with a precision balance



Model	Weighing capacity	Readability	Weighing surface (d)	Price
PCB 200-3	200 g	0.001 g	82 mm	€ 310.-
Suitable Accessories: <ul style="list-style-type: none">• DAkkS calibration certificate 963-127 - 963-127 - €103.-				



<https://partner.kern-sohn.com/Segment/B0CD4>

* DAkkS calibrations are also possible on site
All prices are net, valid until 09.12.2025 plus shipping cost

Personal consultation from our product specialist:



Aleksandar Delic
Productspecialist
Tel : +49 7433 9933-120
aleksandar.delic@kern-sohn.com

Application number: EN - 1434

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com