

RHEUMATOLOGISTS
WEIGH A COTTON BALL BEFORE AND AFTER CHEWING TO DETERMINE THE AMOUNT OF SALIVA

PRECISION BALANCES | KERN 440

To measure the amount of saliva in "Sjögren's syndrome" (reduced tear production and saliva production), a cotton ball is weighed before and after chewing. The measurement is carried out with a precision scale and is an initial check that makes diagnosis easier.



Model	Weighing capacity	Readability	Weighing surface (d)	Price
440-33N	200 g	0.01 g	105 mm	€ 300.-
Suitable Accessories: <ul style="list-style-type: none">• DAkkS calibration certificate 963-127 - 963-127 - €103.-				



<https://partner.kern-sohn.com/Segment/B0FB4>

* DAkkS calibrations are also possible on site
All prices are net, valid until 17.01.2026 plus shipping cost

Personal consultation from our product specialist:



Aleksandar Delic
Productspecialist
Tel : +49 7433 9933-120
aleksandar.delic@kern-sohn.com

Application number: EN - 192 1

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com