

PRODUCTION OF BODY COMPATIBLE SURGICAL MATERIALS
DETERMINATION OF THE GRAIN SIZE OF POWDER PARTICLES BY WEIGHING

PRECISION BALANCES | KERN PNJ

Fast and accurate weighing of the grain size of powder articles in the sieving process. In this way, prescribed limit values for grain sizes can be determined and poor production batches can be avoided. Invalid search result.



Model	Weighing capacity	Readability	Weighing surface (d)	Price
PNJ 3000-2M	3200 g	0.01 g	190×190 mm	€ 890,-
Suitable Accessories: <ul style="list-style-type: none">• DAkkS calibration certificate 963-127 - 963-127 - €103,-				



<https://partner.kern-sohn.com/Segment/S10011>

* DAkkS calibrations are also possible on site
All prices are net, valid until 09.12.2025 plus shipping cost

Personal consultation from our product specialist:



Aleksandar Delic
Productspecialist
Tel : +49 7433 9933-120
aleksandar.delic@kern-sohn.com

Application number: EN - 193

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com