

INTERNISTS
PRECISE WEIGHING OF COTTON BALLS TO DETERMINE SALIVA PRODUCTION (SAXON TEST)

PRECISION BALANCES | KERN 440

To measure saliva production in patients with Sjögren's syndrome, a cotton ball is precisely weighed before and after chewing. A precision scale records the small differences in weight in order to accurately determine the amount of saliva stimulated by chewing.



Model	Weighing capacity	Readability	Weighing surface (d)	Price
440-43N	400 g	0.1 g	130×130 mm	€ 220.-
Suitable Accessories: <ul style="list-style-type: none">• DAkkS calibration certificate 963-127 - 963-127 - €103.-				



<https://partner.kern-sohn.com/Segment/B0FD5>

* DAkkS calibrations are also possible on site
All prices are net, valid until 07.01.2026 plus shipping cost

Personal consultation from our product specialist:



Aleksandar Delic
Productspecialist
Tel : +49 7433 9933-120
aleksandar.delic@kern-sohn.com

Application number: EN - 222 1

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com