

FORENSIC DEPARTMENTS
WEIGHT RECORDING OF AN ORGAN (LIVER, HEART, ETC.) DURING AUTOPSY

PRECISION BALANCES | KERN 440

During an autopsy, the organ that is removed should be weighed directly on the autopsy table, as it could lose weight due to dehydration if left lying around for too long. The reason for recording the organ weight is that there may be a connection between organ weight and the manner of death.



Model	Weighing capacity	Readability	Weighing surface (d)	Price
440-53N	6000 g	1 g	150×170 mm	€ 230.-
Suitable Accessories: <ul style="list-style-type: none">• DAkkS calibration certificate 963-128 - 963-128 - €124.-				



<https://partner.kern-sohn.com/Segment/S10799>

* DAkkS calibrations are also possible on site
All prices are net, valid until 17.01.2026 plus shipping cost

Personal consultation from our product specialist:



Aleksandar Delic
Productspecialist
Tel : +49 7433 9933-120
aleksandar.delic@kern-sohn.com

Application number: EN - 2407

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com