



ENGINEERING OFFICES FOR COATING AND CORROSION PROTECTION WEIGHT CONTROL OF OIL LAYERS ON METAL SAMPLES

ANALYTICAL BALANCES | KERN ABS

In the surface treatment laboratory, metal plates are coated with precisely measured amounts of oil to test corrosion protection and coating properties. An analytical balance allows for precise monitoring of the applied oil layer – this is the only way to ensure reliable and reproducible metal protection.



Model	Weighing capacity	Readability	Price
ABS 320-4N	320 g	0.0001 g	€ 1,550.-
Suitable Accessories: <ul style="list-style-type: none"> • DAkkS calibration certificate 963-101 - 963-101 - €200.- 			



<https://partner.kern-sohn.com/Segment/S10840>

* DAkkS calibrations are also possible on site
All prices are net, valid until 14.11.2025 plus shipping cost

Personal consultation from our product specialist:



Aleksandar Delic
Productspecialist
Tel : +49 7433 9933-120
aleksandar.delic@kern-sohn.com

Application number: EN - 2496

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com