



MANUFACTURE OF STENTS, STENT GRAFTS AND OTHER VASCULAR PROSTHESES, SURGICAL TRANSPLANTS
CHECK WEIGHING OF THE STENT BODY AFTER ELECTROPOLISHING

ANALYTICAL BALANCES | KERN ABT

The analytical balance is essential for inspecting the stent base bodies both after electropolishing and before final assembly. After electropolishing, the mass loss is measured to provide information about the smoothness and roundness of the stents. Before final assembly, the balance enables an intermediate inspection, during which stents that are too large or too small are sorted out to ensure that only the highest quality products are processed.



Model	Weighing capacity	Readability	Price
ABT 120-5DNM	42 g; 120 g	0.00001 g; 0.0001 g	€ 2,550.-
Suitable Accessories: <ul style="list-style-type: none">• DAkkS calibration certificate 963-101 - 963-101 - €200.-			



<https://partner.kern-sohn.com/Segment/S10031>

* DAkkS calibrations are also possible on site
All prices are net, valid until 27.12.2025 plus shipping cost

Personal consultation from our product specialist:



Aleksandar Delic
Productspecialist
Tel : +49 7433 9933-120
aleksandar.delic@kern-sohn.com

Application number: EN - 334

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com

