

INTERNISTS
URINE TEST FOR DOPING CONTROL WITH A REFRACTOMETER

DIGITAL REFRACTOMETERS | KERN ORM-UN

The measurement of the specific urine weight using a refractometer is the validated method in doping control. The refraction of light in the urine is measured, which is directly dependent on the number of particles dissolved in the urine. A very simple, time-saving and elegant method for urine diagnosis



Model	Scales	Scale Range	ATC	Price
ORM 1UN	Refractive index Brix Urine (spec. gravity) Serum proteins	1,3330 nD - 1,4200 nD 0,0 %- 50,0 % 1,000 sgU - 1,050 sgU 0,0 g/dl - 12,0 g/dl		€ 350.-



<https://partner.kern-sohn.com/Segment/B0FD5>

* DAkkS calibrations are also possible on site
All prices are net, valid until 18.01.2026 plus shipping cost

Personal consultation from our product specialist:



Ralf Gutbrod
Productspecialist
Tel : +49 7433 9933-306
ralf.gutbrod@kern-sohn.com

Application number: EN - 873

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com