

DAIRIES

DENSITY TESTING OF WHEY WITH A REFRACTOMETER

DIGITAL REFRACTOMETERS | KERN ORM-R

To ensure the quality of whey powder, the composition of the whey is checked with a refractometer before drying. A Brix value that is too low leads to a dusty end product, so control is crucial for a high-quality end product.



Model	Scales	Scale Range	ATC	Price
ORM 1RS	Refractive index Brix	1,3330 nD - 1,5177 nD 0,0 % - 90,0 %		€ 420,-
Suitable Accessories: <ul style="list-style-type: none"> • Factory calibration certificate 961-290 - 961-290 - €135,- 				



<https://partner.kern-sohn.com/Segment/B00C0>

* DAkkS calibrations are also possible on site
All prices are net, valid until 01.12.2025 plus shipping cost

Personal consultation from our product specialist:



Ralf Gutbrod
Productspecialist
Tel : +49 7433 9933-306
ralf.gutbrod@kern-sohn.com

Application number: EN - 890

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com