

OTOLOGISTS  
WEIGHING A COTTON BALL BEFORE AND AFTER CHEWING TO DIAGNOSE "SJOGREN'S SYNDROME"

PRECISION BALANCES | KERN PCB

To measure the amount of saliva in "Sjögren's Syndrome" (reduced production of tears and saliva), a cotton ball is weighed before and after chewing. The measurement is carried out with a precision scale, which records these small weights very precisely and thus helps to make a diagnosis. Invalid search result.



Model	Weighing capacity	Readability	Weighing surface (d)	Price
PCB 300-3	360 g	0.001 g	82 mm	€ 320.-
<b>Suitable Accessories:</b> <ul style="list-style-type: none"><li>• DAkkS calibration certificate 963-127 - 963-127 - €103.-</li></ul>				



<https://partner.kern-sohn.com/Segment/B0FD2>

\* DAkkS calibrations are also possible on site  
All prices are net, valid until 19.12.2025 plus shipping cost

Personal consultation from our product specialist:



Aleksandar Delic  
Productspecialist  
Tel : +49 7433 9933-120  
aleksandar.delic@kern-sohn.com

Application number: EN - 967

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen, Germany  
info@kern-sohn.com  
[www.kern-sohn.com](http://www.kern-sohn.com)