

OTOLOGISTS

WEIGHING A COTTON BALL BEFORE AND AFTER CHEWING TO DIAGNOSE "SJOGREN'S SYNDROME"

PRECISION BALANCES | KERN PCB

To measure the amount of saliva in "Sjögren's Syndrome" (reduced production of tears and saliva), a cotton ball is weighed before and after chewing. The measurement is carried out with a precision scale, which records these small weights very precisely and thus helps to make a diagnosis. Invalid search result.



Model	Weighing capacity	Readability	Weighing surface (d)	Price
PCB 300-3	360 g	0.001 g	82 mm	€ 320.-
Suitable Accessories:				
<ul style="list-style-type: none"> • DAkkS calibration certificate 963-127 - 963-127 - €103.- 				



<https://partner.kern-sohn.com/Segment/B0FD2-967>

* DAkkS calibrations are also possible on site
All prices are net, valid until 24.10.2025 plus shipping cost

Personal consultation from our product specialist:



Aleksandar Delic
Productspecialist
Tel : +49 7433 9933-120
aleksandar.delic@kern-sohn.com

Application number: EN - 967

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com