



BEEKEEPINGS
WEIGHING HONEY

BENCH SCALES | KERN WTB

This handy table scale with a readability of 0.5 g is ideal for accurately filling honey into 30 g jars. It meets the calibration requirements of food law, thus enabling reliable, legally compliant filling.



Model	Weighing capacity	Readability	Price
WTB 1K-4NM	1500 g	0.5 g	€ 300.-
Suitable Accessories: <ul style="list-style-type: none">• Verification / Conformity assessment according to NAWI (2014/31/EU) 965-227 - 965-227 - €65.-• Rechargeable battery pack WTB-A01N - WTB-A01N - €38.-			



<https://partner.kern-sohn.com/Segment/B006A>

* DAkkS calibrations are also possible on site
All prices plus VAT and shipping costs valid until 02.02.2026

For expert advice, please contact:



Jürgen Fünfer
Productspecialist
Tel : +49 7433 9933-245
juergen.fuenfer@kern-sohn.com

Application number: EN - 1038

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com