



BAKERY FIRMS
WEIGHING OUT SMALL AMOUNTS OF SOURDOUGH TO MAKE A TEST DOUGH

PRECISION BALANCES | KERN 440

The KERN precision balance is specially designed for the accurate measurement of small quantities and is therefore ideal for the development of new baked goods or special recipes.



Model	Weighing capacity	Readability	Weighing surface (d)	Price
440-49N	4000 g	0.1 g	150×170 mm	€ 350,-
Suitable Accessories: <ul style="list-style-type: none">• DAkkS calibration certificate 963-127 - 963-127 - €103,-				



<https://partner.kern-sohn.com/Segment/B00AB>

* DAkkS calibrations are also possible on site
All prices are net, valid until 30.11.2025 plus shipping cost

Personal consultation from our product specialist:



Aleksandar Delic
Productspecialist
Tel : +49 7433 9933-120
aleksandar.delic@kern-sohn.com

Application number: EN - 1046

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com