

BOTTLING AND PACKING COMPANIES
CONTROLLING WEIGHING OF INDIVIDUAL CONTAINERS DURING THE FILLING PROCESS

BENCH SCALES | KERN GAB

Companies use the expertise of a contract bottler to outsource the filling and packaging of their products. The different containers are filled and weighed by hand so that the filling quantity corresponds exactly to the quantity specified on the sticker. A scale with official approval is used for this.



Model	Weighing capacity	Readability	Price
GAB 6K1DNM	3 kg; 6 kg	0.001 kg; 0.002 kg	€ 310.-
Suitable Accessories: <ul style="list-style-type: none">• Verification / Conformity assessment according to NAWI (2014/31/EU) 965-228 - 965-228 - €84.-			



<https://partner.kern-sohn.com/Segment/B129A>

* DAkkS calibrations are also possible on site
All prices are net, valid until 06.01.2026 plus shipping cost

Personal consultation from our product specialist:



Jürgen Fünfer
Productspecialist
Tel : +49 7433 9933-245
juergen.fuenfer@kern-sohn.com

Application number: EN - 1115

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com