


ENGRAVING COMPANY
GOLD SCALE FOR WEIGHING JEWELLERY

PRECISION BALANCES | KERN PBJ

This scale with 6 kg capacity and 0.01g readability is used to weigh gold jewellery. The scale can be calibrated and is therefore suitable for commercial use. Robust and stable thanks to the metal housing



Model	Weighing capacity	Readability	Weighing surface (d)	Price
PBJ 6200-2M	6200 g	0.01 g	180×190 mm	€ 2,300.-
Suitable Accessories: <ul style="list-style-type: none">• Verification / Conformity assessment according to NAWI (2014/31/EU) 965-202 – 965-202 – €150.-				



<https://partner.kern-sohn.com/Segment/B05DF>

* DAkkS calibrations are also possible on site
All prices are net, valid until 13.11.2025 plus shipping cost

Personal consultation from our product specialist:



Aleksandar Delic
Productspecialist
Tel : +49 7433 9933-120
aleksandar.delic@kern-sohn.com

Application number: EN - 1153

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com