


CREMATORIUM
WEIGHING OF COFFINS

FLOOR SCALES | KERN BIC

Upon arrival at the crematorium, the coffin is weighed on floor scales. The weight decides in which cremation facility the cremation takes place. Thus, the safety of the cremation is guaranteed at all times



Model	Weighing capacity	Readability	Price
BIC 600K-1	300 kg; 600 kg	0.1 kg; 0.2 kg	€ 1,340.-
Suitable Accessories: <ul style="list-style-type: none">● Pit frame BIC-A05 - BIC-A05 - €460.-● Pair of base plates BIC-A07 - BIC-A07 - €44.-● Access ramp BIC-A02, 1200×1000×108 mm - BIC-A02 - €530.-			

* DAkkS calibrations are also possible on site
All prices are net, valid until 06.11.2025 plus shipping cost



<https://partner.kern-sohn.com/Segment/S10536>

Personal consultation from our product specialist:



Aleksandar Delic
Productspecialist
Tel : +49 7433 9933-120
aleksandar.delic@kern-sohn.com

Application number: EN - 1286

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com