



HIGH-PERFORMANCE CERAMICS MANUFACTURER
DETERMINE THE MOISTURE CONTENT OF RAW MATERIALS

MOISTURE ANALYSER | KERN DBS

Precise determination of the moisture content of the mixture of clay, loam and water (so-called slip). The correct moisture content of the raw materials used is a prerequisite for obtaining a flawless product



Model	Weighing capacity	Readability	Price
DBS 60-3	60 g	0.001 g	€ 1,650.-
Suitable Accessories: <ul style="list-style-type: none">• Testing individual substances (Moisture Analysers) 965-243 - 965-243 - €105.-• Sample plates aluminium MLB-A01A - MLB-A01A - €45.-			



<https://partner.kern-sohn.com/Segment/S10431>

* DAkkS calibrations are also possible on site
All prices are net, valid until 06.01.2026 plus shipping cost

Personal consultation from our product specialist:



Aleksandar Delic
Productspecialist
Tel : +49 7433 9933-120
aleksandar.delic@kern-sohn.com

Application number: EN - 136 1

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com