



JEWELRY MANUFACTURER
CHECKING THE WEIGHT OF A FINISHED JEWELRY RING

PRECISION BALANCES | KERN PNJ

Jewelery making is solid craftsmanship and before a ring leaves the manufactory there is a strict quality control. Here the inspector checks whether each stone is firmly set, the quality is impeccable and the weight of the ring corresponds to the specification. A precision scale is used to check the weight, which is ideally suited for these small measurements



Model	Weighing capacity	Readability	Weighing surface (d)	Price
PNJ 3000-2M	3200 g	0.01 g	190×190 mm	€ 890.-
Suitable Accessories: <ul style="list-style-type: none">• Verification / Conformity assessment according to NAWI (2014/31/EU) 965-216 - 965-216 - €85.-				



<https://partner.kern-sohn.com/Segment/B0A30>

* DAkkS calibrations are also possible on site
All prices are net, valid until 10.11.2025 plus shipping cost

Personal consultation from our product specialist:



Aleksandar Delic
Productspecialist
Tel : +49 7433 9933-120
aleksandar.delic@kern-sohn.com

Application number: EN - 1363

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com