

MANUFACTURERS OF TINS OF MEAT AND SAUSAGE
WEIGHING COMPONENTS FOR A SPICE MIX

BENCH SCALES | KERN WTB

Butchers use salt, pepper, allspice, ginger, cardamom, garlic, herbs and much more to season different types of canned sausage. The spices have to be weighed very carefully, as the sausage must always have the same standard of taste in order to satisfy the consumer.



Model	Weighing capacity	Readability	Price
WTB 1K-4N	1500 g	0.2 g	€ 285.-
Suitable Accessories: <ul style="list-style-type: none">DAkkS calibration certificate 963-127 - 963-127 - €103.-			



<https://partner.kern-sohn.com/Segment/B00A7>

* DAkkS calibrations are also possible on site
All prices are net, valid until 15.12.2025 plus shipping cost

Personal consultation from our product specialist:



Jürgen Fünfer
Productspecialist
Tel : +49 7433 9933-245
juergen.fuenfer@kern-sohn.com

Application number: EN - 1372

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com