

REFRIGERATED COUNTERS READY MEALS MANUFACTURER
WEIGH THE INDIVIDUAL COMPONENTS (VEGETABLES, SALADS, MEAT, ETC.) FOR A READY MEAL

BENCH SCALES | KERN FOB

When preparing the ultra-fresh ready meals (chilled food), all components such as vegetables, salads, meat or fish are precisely weighed with scales. Accurate portioning is important for customer satisfaction.



Model	Weighing capacity	Readability	Price
FOB 3K1	3000 g	1 g	€ 240.-
Suitable Accessories: <ul style="list-style-type: none">• DAkkS calibration certificate 963-127 - 963-127 - €103.-			



<https://partner.kern-sohn.com/Segment/S10545>

* DAkkS calibrations are also possible on site
All prices are net, valid until 14.12.2025 plus shipping cost

Personal consultation from our product specialist:



Jürgen Fünfer
Productspecialist
Tel : +49 7433 9933-245
juergen.fuenfer@kern-sohn.com

Application number: EN - 1386

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com