



PARQUET FLOORERS  
FRICTION TEST ON FRESHLY LAID FLOORS TO TEST ANTI-SLIP PRODUCTS

DIGITAL FORCE-MEASURING DEVICES (EXTERNAL SENSOR) | KERN FL

The effectiveness of anti-slip products can be checked with a friction test. A plate is pulled over a surface and the pulling force that needs to be applied is checked before and after the anti-slip treatment. A very fast and efficient testing method to get results. Possible risks due to product liability can be reduced in this way



| Model   | Max. Force | Readability | Price     |
|---|------------|-------------|-----------|
| FL 20K  | 20 kN      | 10 N        | € 1,290.- |
| <b>Suitable Accessories:</b> <ul style="list-style-type: none"><li>• Factory calibration certificate 961-164 (tension   N) - 961-164 - €410.-</li></ul> |            |             |           |



<https://partner.kern-sohn.com/Segment/B096D-1530>

\* DAkkS calibrations are also possible on site  
All prices are net, valid until 25.10.2025 plus shipping cost

Personal consultation from our product specialist:



Irmi Russo Productspecialist  
Tel : +49 7433 9933-208  
russo@kern-sohn.com

Application number: EN - 1530  
**KERN & SOHN GmbH**  
Ziegelei 1  
72336 Balingen, Germany  
info@kern-sohn.com  
[www.kern-sohn.com](http://www.kern-sohn.com)