



PRECIOUS AND SEMI-PRECIOUS STONES EXPORTERS
CHECKING THE QUALITY OF GEMSTONES USING THE REFRACTIVE INDEX

ANALOGUE REFRACTOMETERS | KERN ORA-G

When exporting gemstones, quality control is an important control process to verify the authenticity of the stone. Each gemstone has its typical refractive index due to its composition and crystalline structure. The refractive index quality control test is a quick and elegant method of inspecting and analyzing the stone



Model	Scales	Scale Range	ATC	Price
ORA 1GG	Refractive index	1,30 nD – 1,81 nD		€ 250.-



<https://partner.kern-sohn.com/Segment/B0A17>

* DAkkS calibrations are also possible on site
All prices plus VAT and shipping costs valid until 28.03.2026

We would be happy to advise you personally:



Ralf Gutbrod
Productspecialist
Tel : +49 7433 9933-306
ralf.gutbrod@kern-sohn.com

Application number: EN - 1532

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com