



PARCEL SERVICES, LOGISTIC, RELOADING DEPOTS
AUTOMATIC WEIGHING OF ALL PACKAGES TO DETERMINE THE SHIPPING WEIGHT

PLATFORM SCALES | KERN IOC

Automatic and accurate weighing of all parcels to determine the shipping weight and to correctly pass on the transport price to the customer. A simple, precise process that is extremely time-saving



Model	Weighing capacity	Readability	Price
IOC 100K-2LM	60 kg; 150 kg	0.02 kg; 0.05 kg	€ 700.-
Suitable Accessories: <ul style="list-style-type: none">• DAkkS calibration certificate 963-129 - 963-129 - €153.-			



<https://partner.kern-sohn.com/Segment/S10020>

* DAkkS calibrations are also possible on site
All prices are net, valid until 11.11.2025 plus shipping cost

Personal consultation from our product specialist:



Jürgen Fünfer
Productspecialist
Tel : +49 7433 9933-245
juergen.fuenfer@kern-sohn.com

Application number: EN - 157

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com