



BEEKEEPER CLUBS
ARTIFICIAL INSEMINATION OF QUEEN BEES UNDER THE MICROSCOPE

STEREO MICROSCOPES | KERN OZM-5

Beekeepers can inseminate queen bees in a targeted manner - via so-called artificial instrumental insemination - in order to achieve certain breeding goals. A microscope is used to hold the bee and insert the sperm. It enables precise, faster and animal-friendly work



Model	Optical System	Tube	Eyepiece	Price
OZM 542	Greenough	Binocular	Eyepiece HSWF 10 x / Ø 23mm with anti-fungus, high eye point	€ 1,390.-



<https://partner.kern-sohn.com/Segment/B1502>

* DAkkS calibrations are also possible on site
 All prices are net, valid until 20.11.2025 plus shipping cost

Personal consultation from our product specialist:

Application number: EN - 1587



Ralf Gutbrod
 Productspecialist
 Tel : +49 7433 9933-306
ralf.gutbrod@kern-sohn.com

KERN & SOHN GmbH
 Ziegelei 1
 72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com