

BEEKEEPING CLUBS  
INSEMINATION OF A QUEEN BEE UNDER THE MICROSCOPE

STEREO MICROSCOPES | KERN OSF-4G

For the artificial insemination of a queen bee, equipment consisting of the actual insemination device for sperm transfer and a stereo microscope for optical control is required,The microscope used for this can also be used to mark a queen bee with a color code,All manipulations are possible under the microscope,A magnification of 10x is usual



Model	Optical System	Tube	Eyepiece	Price
OSF 438	Greenough	Binocular	Eyepiece WF 10x / Ø 20mm with anti-fungus	€ 390.-



<https://partner.kern-sohn.com/Segment/B15DF>

\* DAkkS calibrations are also possible on site  
All prices are net, valid until 06.01.2026 plus shipping cost

Personal consultation from our product specialist:



**Ralf Gutbrod**  
Productspecialist  
Tel : +49 7433 9933-306  
ralf.gutbrod@kern-sohn.com

Application number: EN - 1623

**KERN & SOHN GmbH**  
Ziegelei 1  
72336 Balingen, Germany  
info@kern-sohn.com  
www.kern-sohn.com