

TROUT FARMS  
DETERMINING THE WEIGHT OF TROUT SOLD DIRECTLY AT THE BREEDING TANK

BENCH SCALES | KERN GAB

Trout are fished directly from the water and weighed to determine the price based on the weight, This type of direct marketing offers the advantage that the customer can get a fish according to his idea (weight).



Model	Weighing capacity	Readability	Price
GAB 15K2DNM	6000 g; 15000 g	2 g; 5 g	€ 310.-
<b>Suitable Accessories:</b> <ul style="list-style-type: none"><li>• Verification / Conformity assessment according to NAWI (2014/31/EU) 965-228 - 965-228 - €84.-</li></ul>			



<https://partner.kern-sohn.com/Segment/B0045>

\* DAkkS calibrations are also possible on site  
All prices are net, valid until 21.11.2025 plus shipping cost

Personal consultation from our product specialist:



Jürgen Fünfer  
Productspecialist  
Tel : +49 7433 9933-245  
[juergen.fuenfer@kern-sohn.com](mailto:juergen.fuenfer@kern-sohn.com)

Application number: EN - 1633

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen, Germany  
[info@kern-sohn.com](mailto:info@kern-sohn.com)  
[www.kern-sohn.com](http://www.kern-sohn.com)