


TEXTILES FOR INTERIOR DESIGN
WEIGHING A PACKAGE READY FOR SHIPMENT TO DETERMINE THE FREIGHT CHARGE

PLATFORM SCALES | KERN IOC

After production of painted, flocked or printed room textiles, they are folded and packaged for shipping. The package is weighed to determine the freight fee and for registration with the parcel service provider. The package weight, which is correctly determined in advance, enables the shipping process to be processed quickly, thereby saving time, avoiding any additional franking and thereby saving costs



Model	Weighing capacity	Readability	Price
IOC 60K-2LM	30 kg; 60 kg	0.01 kg; 0.02 kg	€ 500.-
Suitable Accessories: <ul style="list-style-type: none">• Verification / Conformity assessment according to NAWI (2014/31/EU) 965-229 - 965-229 - €110.-			



<https://partner.kern-sohn.com/Segment/B023F>

* DAkkS calibrations are also possible on site
All prices are net, valid until 22.11.2025 plus shipping cost

Personal consultation from our product specialist:



Jürgen Fünfer
Productspecialist
Tel : +49 7433 9933-245
juergen.fuenfer@kern-sohn.com

Application number: EN - 1790

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com