

INDOOR FARMING
CONTROLLING WEIGHING OF BERRY PACKAGES USING A CALIBRATED SCALE

BENCH SCALES | KERN FFN

After harvest, berries grown in indoor farms are packaged in 250g or 500g containers. The weight is carefully controlled as filling too high would reduce the return. Insufficient filling would result in complaints.



Model	Weighing capacity	Readability	Price
FFN 6K2IPM	6 kg	0.002 kg	€ 340.-
Suitable Accessories: <ul style="list-style-type: none">• Verification / Conformity assessment according to NAWI (2014/31/EU) 965-228 - 965-228 - €84.-			



<https://partner.kern-sohn.com/Segment/S10460>

* DAkkS calibrations are also possible on site
All prices are net, valid until 17.11.2025 plus shipping cost

Personal consultation from our product specialist:



Jürgen Fünfer
Productspecialist
Tel : +49 7433 9933-245
juergen.fuenfer@kern-sohn.com

Application number: EN - 1800

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com