



GEMSTONE CUTTER

QUALITY TESTING OF A GEMSTONE BY THE REFRACTIVE INDEX

ANALOGUE REFRACTOMETERS | KERN ORA-G

Gemstone quality testing is an important control process to verify the authenticity of a stone. Quality control testing using a refractometer, which determines the refractive index, is a quick and elegant method for assessing and analyzing the stone and ultimately determining its value.



Model	Scales	Scale Range	ATC	Price
ORA 1GG	Refractive index	1,30 nD – 1,81 nD		€ 250.-



<https://partner.kern-sohn.com/Segment/S10281>

* DAkkS calibrations are also possible on site
All prices plus VAT and shipping costs valid until 28.03.2026

We would be happy to advise you personally:



Ralf Gutbrod
Productspecialist
Tel : +49 7433 9933-306
ralf.gutbrod@kern-sohn.com

Application number: EN - 197 1

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com