

## FLOOR COVERS MANUFACTURER TESTING THE TEAR RESISTANCE OF A FLOOR COVERING

### TEST STANDS (MOTORISED) | KERN TVO

A plastic floor made of granules and composite materials etc. is subjected to various laboratory tests after production, either randomly or depending on the batch. The tear resistance is also tested. This test is an important guarantee for the safety of the floor covering and avoids subsequent damage or complaints.



Model	Max. Force	Price
TVO 500N300	500 N	€ 2,500.-
<b>Suitable Accessories:</b> <ul style="list-style-type: none"><li>• Max. Force 500 N, Readability 0.2 N - FL 500 - €570.-</li><li>• Factory calibration certificate 961-361 (tension, compression   N) - 961-361 - €340.-</li></ul>		



<https://partner.kern-sohn.com/Segment/B0BE6>

\* DAkkS calibrations are also possible on site  
All prices are net, valid until 07.01.2026 plus shipping cost

Personal consultation from our product specialist:



**Irmi Russo** Productspecialist  
Tel : +49 7433 9933-208  
[russo@kern-sohn.com](mailto:russo@kern-sohn.com)

Application number: EN - 2047

**KERN & SOHN GmbH**  
Ziegelei 1  
72336 Balingen, Germany  
[info@kern-sohn.com](mailto:info@kern-sohn.com)  
[www.kern-sohn.com](http://www.kern-sohn.com)

