

SEED BANKS

PRECISE WEIGHING OF SEEDS FOR SHIPPING TO CUSTOMERS

PRECISION BALANCES | KERN PCJ

In seed breeding companies, dried seeds are precisely weighed after harvest in order to send them to hobby gardeners and specialist companies. Precision scales are used that can accurately record even the smallest weight differences of 0.001g to 50g. This prevents inaccuracies that could lead to economic disadvantages for both the breeder and the customer.



Model	Weighing capacity	Readability	Weighing surface (d)	Price
PCJ 600-2M	600 g	0.01 g	131×131 mm	€ 420.-
Suitable Accessories:				
<ul style="list-style-type: none"> • Verification / Conformity assessment according to NAWI (2014/31/EU) 965-216 - 965-216 - €85.- 				



<https://partner.kern-sohn.com/Segment/B0020-2117>

* DAkkS calibrations are also possible on site
All prices are net, valid until 03.10.2025 plus shipping cost

Personal consultation from our product specialist:

Application number: EN - 2117



Aleksandar Delic
Productspecialist
Tel : +49 7433 9933-120
aleksandar.delic@kern-sohn.com

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com