

SEED BANKS  
PRECISE WEIGHING OF SEEDS FOR SHIPPING TO CUSTOMERS

PRECISION BALANCES | KERN PCJ

In seed breeding companies, dried seeds are precisely weighed after harvest in order to send them to hobby gardeners and specialist companies. Precision scales are used that can accurately record even the smallest weight differences of 0.001g to 50g. This prevents inaccuracies that could lead to economic disadvantages for both the breeder and the customer.



Model	Weighing capacity	Readability	Weighing surface (d)	Price
PCJ 600-2M	600 g	0.01 g	131×131 mm	€ 420.-
<b>Suitable Accessories:</b> <ul style="list-style-type: none"><li>• Verification / Conformity assessment according to NAWI (2014/31/EU) 965-216 - 965-216 - €85.-</li></ul>				



<https://partner.kern-sohn.com/Segment/B0020>

\* DAkkS calibrations are also possible on site  
All prices are net, valid until 30.11.2025 plus shipping cost

Personal consultation from our product specialist:



Aleksandar Delic  
Productspecialist  
Tel : +49 7433 9933-120  
aleksandar.delic@kern-sohn.com

Application number: EN - 2 1 1 7

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen, Germany  
info@kern-sohn.com  
www.kern-sohn.com