



MUSHROOMS WHOLESALE
WEIGHING MUSHROOMS DURING HARVEST

BENCH SCALES | KERN ECB

Since mushrooms are very sensitive, they are only picked once during harvesting and immediately placed in the sales bowl for a so-called pre-test of the weight. Packages of 250g to 3kg are weighed. Careful handling during harvesting protects the mushroom body, avoids bruising and possible rotting during transport.



Model	Weighing capacity	Readability	Price
ECE 10K-3N	10 kg	0.005 kg	€ 120.-
Suitable Accessories: <ul style="list-style-type: none">• DAkkS calibration certificate 963-128 - 963-128 - €124.-			



<https://partner.kern-sohn.com/Segment/B01BA>

* DAkkS calibrations are also possible on site
All prices plus VAT and shipping costs valid until 26.01.2026

For expert advice, please contact:



Jürgen Fünfer
Productspecialist
Tel : +49 7433 9933-245
juergen.fuenfer@kern-sohn.com

Application number: EN - 2 147

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com