



**SHRIMP FARMS**  
PRECISE WEIGHING OF FEED IN SHRIMP FARMING

PRECISION BALANCES | KERN PCJ

In shrimp farms, precise dosing of protein sources such as freeze-dried bloodworms is essential. A precision scale ensures that the amount of feed is precisely adjusted to the age and number of animals. This contributes to the optimal health and growth of the shrimp.



Model	Weighing capacity	Readability	Weighing surface (d)	Price
PCJ 600-2M	600 g	0.01 g	131×131 mm	€ 420.-



<https://partner.kern-sohn.com/Segment/S10809-2165>

\* DAkkS calibrations are also possible on site  
All prices are net, valid until 14.09.2025 plus shipping cost

Personal consultation from our product specialist:

Application number: EN - 2165



**Aleksandar Delic**  
Productspecialist  
Tel : +49 7433 9933-120  
aleksandar.delic@kern-sohn.com

**KERN & SOHN GmbH**  
Ziegelei 1  
72336 Balingen, Germany  
info@kern-sohn.com  
www.kern-sohn.com