



SHRIMP FARMS
PRECISE WEIGHING OF FEED IN SHRIMP FARMING

PRECISION BALANCES | KERN PCJ

In shrimp farms, precise dosing of protein sources such as freeze-dried bloodworms is essential. A precision scale ensures that the amount of feed is precisely adjusted to the age and number of animals. This contributes to the optimal health and growth of the shrimp.



Model	Weighing capacity	Readability	Weighing surface (d)	Price
PCJ 600-2M	600 g	0.01 g	131×131 mm	€ 420.-



<https://partner.kern-sohn.com/Segment/S10809-2165>

* DAkkS calibrations are also possible on site
All prices are net, valid until 14.09.2025 plus shipping cost

Personal consultation from our product specialist:

Application number: EN - 2165



Aleksandar Delic
Productspecialist
Tel : +49 7433 9933-120
aleksandar.delic@kern-sohn.com

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com