PROFESSIONAL MEASURING





SHRIMP FARMS PRECISE WEIGHING OF FEED IN SHRIMP FARMING

PRECISION BALANCES | KERN PCJ

In shrimp farms, precise dosing of protein sources such as freeze-dried bloodworms is essential. A precision scale ensures that the amount of feed is precisely adjusted to the age and number of animals. This contributes to the optimal health and growth of the shrimp.



Model	Weighing capacity	Readability	Weighing surface (d)	Price
PCJ 600-2M	600 g	0.01 g	131×131 mm	€ 420



https://partner.kern-sohn.com/Segment/S10809

* DAkkS calibrations are also possible on site All prices plus VAT and shipping costs valid until 30.04.2026

We would be happy to advise you personally:



Aleksandar Delic
Productspecialist
Tel: +49 7433 9933-120
aleksandar.delic@kern-sohn.com

Application number: EN - 2165

KERN & SOHN GmbH

Tiogoloi 1

Ziegelei 1 72336 Balingen, Germany info@kern-sohn.com www.kern-sohn.com