



MEDICAL INSTRUMENTS MANUFACTURER
INVENTORY RECORDING OF ASSEMBLY PARTS IN THE WAREHOUSE WITH A COUNTING SCALE

COUNTING SCALES | KERN CDS

The counting scale enables quick and reliable recording of small parts in the warehouse, such as for the assembly of medical specula. Using the counting scale makes inventory recording more efficient by avoiding errors in manual counting and significantly speeding up the process. This saves both time and money.



Model	Weighing capacity	Readability	Counting resolution	Price
CDS 30K0.1	30 kg	0.0001 kg	300000	€ 670.-
Suitable Accessories:				
<ul style="list-style-type: none"> • DAkkS calibration certificate 963-128 - 963-128 - €118.- 				



<https://partner.kern-sohn.com/Segment/B0905-2183>

* DAkkS calibrations are also possible on site
 All prices are net, valid until 03.08.2025 plus shipping cost

Personal consultation from our product specialist:



Jürgen Fünfer
 Productspecialist
 Tel : +49 7433 9933-245
juergen.fuenfer@kern-sohn.com

Application number: EN - 2183

KERN & SOHN GmbH
 Ziegelei 1
 72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com