

MANUFACTURERS OF ANIMAL INSECT FEED
CONTROLLED WEIGHING OF READY-PACKED CONTAINERS WITH ANIMAL FEED

BENCH SCALES | KERN WTB

Once the feed made from animal meal or insect meal has been mixed, containers of 50g-2kg are packaged and weighed to check. A table scale is used for this, which should be calibrated. This check prevents too little being packaged, which could lead to complaints. However, if too much is packaged, this reduces the return.



Model	Weighing capacity	Readability	Price
WTB 3K-3NM	3000 g	1 g	€ 300.-
Suitable Accessories: <ul style="list-style-type: none">• DAkkS calibration certificate 963-127 - 963-127 - €103.-			



<https://partner.kern-sohn.com/Segment/B0097>

* DAkkS calibrations are also possible on site
All prices are net, valid until 21.11.2025 plus shipping cost

Personal consultation from our product specialist:



Jürgen Fünfer
Productspecialist
Tel : +49 7433 9933-245
juergen.fuenfer@kern-sohn.com

Application number: EN - 2363

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com