

LADIES UNDERWEAR
WEIGHING DYE AND BINDING CARRIER FOR A DYE BATH

BENCH SCALES | KERN FCB

For the textile dye bath, two components—dye and binding agent—are weighed. Accurate weighing is crucial to ensure that individual batches do not exhibit color variations. As a producer, this allows you to meet your customers' needs and expectations and avoid complaints.



Model	Weighing capacity	Readability	Price
FCB 12K1	12 kg	0.001 kg	€ 220.-
Suitable Accessories: <ul style="list-style-type: none">• DAkkS calibration certificate 963-128 - 963-128 - €124.-			



<https://partner.kern-sohn.com/Segment/B028C-2433>

* DAkkS calibrations are also possible on site
All prices are net, valid until 24.10.2025 plus shipping cost

Personal consultation from our product specialist:



Jürgen Fünfer
Productspecialist
Tel : +49 7433 9933-245
juergen.fuenfer@kern-sohn.com

Application number: EN - 2433

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com