

DISINFECTANTS MANUFACTURERS
TESTING THE WALL THICKNESS OF A DISINFECTANT CONTAINER

MATERIAL THICKNESS GAUGES | KERN TO

Before filling a disinfectant container, the material thickness should be checked to ensure the container can withstand the chemical requirements. A material thickness gauge can be used to precisely determine the wall thickness—with or without a coating. This safety check prevents material failure and ensures reliable storage!



Model	Max. Pulse-Echo	Price
TO 100-0.01EE	600 mm	€ 1,490.-
Suitable Accessories: <ul style="list-style-type: none">• Factory calibration certificate 961-113 - 961-113 - €167.-		



<https://partner.kern-sohn.com/Segment/B03E9>

* DAkkS calibrations are also possible on site
All prices are net, valid until 06.11.2025 plus shipping cost

Personal consultation from our product specialist:



Irmi Russo Productspecialist
Tel : +49 7433 9933-208
russo@kern-sohn.com

Application number: EN - 244 1

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com