

LADIES NIGHTWEAR
WEIGHING DYE AND BINDER TO MIX A DYE BATH FOR TEXTILES

PLATFORM SCALES | KERN SFE

For the textile dye bath, two components—dye and binding agent—are weighed. Accurate weighing is crucial to ensure that individual batches do not exhibit color variations. As a producer, this allows you to meet your customers' needs and expectations and avoid complaints.



Model	Weighing capacity	Readability	Price
SFE 100K-2XLNM	150 kg	0.05 kg	€ 930.-
Suitable Accessories: <ul style="list-style-type: none">• DAkkS calibration certificate 963-129 - 963-129 - €153.-			



<https://partner.kern-sohn.com/Segment/B028E>

* DAkkS calibrations are also possible on site
All prices are net, valid until 10.11.2025 plus shipping cost

Personal consultation from our product specialist:



Jürgen Fünfer
Productspecialist
Tel : +49 7433 9933-245
juergen.fuenfer@kern-sohn.com

Application number: EN - 25 10

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com