

LADIES NIGHTWEAR
WEIGHING DYE AND BINDER TO MIX A DYE BATH FOR TEXTILES

PLATFORM SCALES | KERN SFE

For the textile dye bath, two components—dye and binding agent—are weighed. Accurate weighing is crucial to ensure that individual batches do not exhibit color variations. As a producer, this allows you to meet your customers' needs and expectations and avoid complaints.



| Model | Weighing capacity | Readability | Price |
|---|-------------------|-------------|---------|
| SFE 100K-2XLNM | 150 kg | 0.05 kg | € 860.- |
| Suitable Accessories: <ul style="list-style-type: none">• DAkkS calibration certificate 963-129 - 963-129 - €161.- | | | |



<https://partner.kern-sohn.com/Segment/B028E>

* DAkkS calibrations are also possible on site
All prices plus VAT and shipping costs valid until 16.07.2026

We would be happy to advise you personally:



Jürgen Fünfer
Productspecialist
Tel : +49 7433 9933-245
juergen.fuenfer@kern-sohn.com

Application number: EN - 25 10

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com