

GOAT CHEESE MANUFACTURER  
PRECISE MEASURING OF RENNET POWDER FOR GOAT CHEESE PRODUCTION

PRECISION BALANCES | KERN 440

To make goat cheese, rennet powder is precisely measured using precision scales. Even the smallest deviations in the amount can affect the consistency, ripeness, and quality of the cheese – therefore, accurate weighing is essential for a consistently high-quality end product!



Model	Weighing capacity	Readability	Weighing surface (d)	Price
440-47N	2000 g	0.1 g	130×130 mm	€ 300.-
<b>Suitable Accessories:</b> <ul style="list-style-type: none"><li>• DAkkS calibration certificate 963-127 - 963-127 - €103.-</li></ul>				



<https://partner.kern-sohn.com/Segment/S10382>

\* DAkkS calibrations are also possible on site  
All prices are net, valid until 09.12.2025 plus shipping cost

Personal consultation from our product specialist:



**Aleksandar Delic**  
Productspecialist  
Tel : +49 7433 9933-120  
aleksandar.delic@kern-sohn.com

Application number: EN - 25 17

**KERN & SOHN GmbH**  
Ziegelei 1  
72336 Balingen, Germany  
info@kern-sohn.com  
[www.kern-sohn.com](http://www.kern-sohn.com)