



LOGISTICS SERVICES
PICKING OF GOODS IN BITO BOXES BY WEIGHING

PLATFORM SCALES | KERN IOC

Warehouse picking of small items from the incoming goods pallet in BITO boxes for storage in the shuttle warehouse (semi-automatic warehouse). Weight test against overload in the box: 30 kg. High profitability by reducing operating costs



Model	Weighing capacity	Readability	Price
IOC 60K-3L	30 kg; 60 kg	0.001 kg; 0.002 kg	€ 510.-



<https://partner.kern-sohn.com/Segment/S10016>

* DAkkS calibrations are also possible on site
All prices plus VAT and shipping costs valid until 26.01.2026

For expert advice, please contact:



Jürgen Fünfer
Productspecialist
Tel : +49 7433 9933-245
juergen.fuenfer@kern-sohn.com

Application number: EN - 307

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com