



LOGISTIC SERVICES
CHECKING THE COMPLETENESS OF THE FILLED BOXES DURING ORDER PICKING BY CHECKWEIGHING

PLATFORM SCALES | KERN IOC

Completeness check of filled boxes during order picking by means of checkweighing. Comparison of the actual weight with the target weight according to the order in the ERP system. Time-saving checking of the boxes ready for dispatch.



Model	Weighing capacity	Readability	Price
IOC 60K-3L	30 kg; 60 kg	0.001 kg; 0.002 kg	€ 510.-



<https://partner.kern-sohn.com/Segment/S10016>

* DAkkS calibrations are also possible on site
All prices are net, valid until 04.11.2025 plus shipping cost

Personal consultation from our product specialist:



Jürgen Fünfer
Productspecialist
Tel : +49 7433 9933-245
juergen.fuenfer@kern-sohn.com

Application number: EN - 308

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com