



MANUFACTURE OF STENTS, STENT GRAFTS AND OTHER VASCULAR PROSTHESES, SURGICAL TRANSPLANTS

TENSILE TESTING OF THE EPTFE VASCULAR PROSTHESES

TEST STANDS (MOTORISED) | KERN TVO

The motorized test rig is used for tensile testing of ePTFE vascular grafts to verify their tear resistance. These tests ensure that the grafts can withstand the high mechanical demands. They safeguard product quality and reduce the risk of product defects.



Model	Max. Force	Price
TVO 2000N500S	2 kN	€ 5,800.-



<https://partner.kern-sohn.com/Segment/S10031>

* DAkkS calibrations are also possible on site
All prices are net, valid until 17.01.2026 plus shipping cost

Personal consultation from our product specialist:



Irmi Russo Productspecialist

Tel : +49 7433 9933-208

russo@kern-sohn.com

Application number: EN - 329

KERN & SOHN GmbH

Ziegelei 1

72336 Balingen, Germany

info@kern-sohn.com

www.kern-sohn.com