



MANUFACTURE OF STENTS, STENT GRAFTS AND OTHER VASCULAR PROSTHESES, SURGICAL TRANSPLANTS
DENSITY DETERMINATION OF THE EPTFE PROSTHESIS MATERIAL

ANALYTICAL BALANCES | KERN ABT

Density determination of the ePTFE prosthesis material in order to guarantee a perfect product during manufacture



Model	Weighing capacity	Readability	Price
ABT 320-4NM	320 g	0.0001 g	€ 2,600.-
Suitable Accessories: <ul style="list-style-type: none">• Density determination set YDB-03 - YDB-03 - €380.-• DAkkS calibration certificate 963-101 - 963-101 - €200.-			



<https://partner.kern-sohn.com/Segment/S10031>

* DAkkS calibrations are also possible on site
All prices are net, valid until 06.12.2025 plus shipping cost

Personal consultation from our product specialist:



Aleksandar Delic
Productspecialist
Tel : +49 7433 9933-120
aleksandar.delic@kern-sohn.com

Application number: EN - 332

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com