


CONTRACT TURNING SHOPS FOR CUSTOMER-SPECIFIC PARTS

WEIGHING OF THE LOADED PALLETS IN THE OUTGOING GOODS AREA

COUNTING SYSTEMS | KERN CCS

Weighing of the loaded pallets in the outgoing goods area. Since empty pallets do not have a standardized dead weight, it is advisable to weigh them individually after loading. This avoids incorrect weight information from the forwarder and possible additional costs



Model	Counting resolution	Readout weight reference scale	Max weight quantity scale	Price
CCS 600K-1S	6000000	100 mg	600 kg	€ 1,840.-
Suitable Accessories: <ul style="list-style-type: none">• DAkkS calibration certificate 962-130-127 - 962-130-127 - €320.-				



<https://partner.kern-sohn.com/Segment/S10003>

* DAkkS calibrations are also possible on site
All prices are net, valid until 07.12.2025 plus shipping cost

Personal consultation from our product specialist:



Jürgen Fünfer
Productspecialist
Tel : +49 7433 9933-245
juergen.fuenfer@kern-sohn.com

Application number: EN - 345

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com