



DEPOSIT THICKNESS GAUGE MANUFACTURER (TESTSEGMENT)
TEST BENCH FOR MULTI-PURPOSE TESTS

TEST STANDS (MOTORISED) | KERN TVM

Super test bench for great multi-purpose tests



| Model | Max. Force | Price |
|---------------|------------|-----------|
| TVM 5000N230N | 5 kN | € 2,740.- |



<https://partner.kern-sohn.com/Segment/B08E7>

* DAkkS calibrations are also possible on site
All prices are net, valid until 09.11.2025 plus shipping cost

Personal consultation from our product specialist:



Irmi Russo Productspecialist
Tel : +49 7433 9933-208
russo@kern-sohn.com

Application number: EN - 384

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com